

Exhibit B



About Guardant360

Guardant360 is a comprehensive liquid profiling test that evaluates 740 genes and includes TMB, MSI status, CHIP filtering, and methylation-based features.

Test Specifications

| Sample input | Sample specifications | Turnaround time |
|--------------------------------|--|--|
| Two 10 mL tubes of whole blood | Ship same or next day at room temperature - do not freeze or refrigerate | 7 days from sample receipt to results* |

Performance Specifications

| Alterations | Limit of Detection at 95% Sensitivity [†] | Specificity | Threshold for Positivity [‡] |
|------------------------|--|-------------|---------------------------------------|
| SNVs | 0.20% | ≥99.9% | ≥0.001% |
| Indels | 0.26% | ≥99.9% | ≥0.01% |
| CNAs | 2.46 copies | ≥99.9% | ≥2.16 copies |
| Fusions/Rearrangements | 0.15% | ≥99.9% | ≥2 unique molecules |
| MSI-High | 0.05% | ≥99.9% | — |
| Tumor Fraction | 0.05% | — | — |
| TMB | ≥0.3% [§] | — | — |
| Promoter Methylation | 1.6% | — | ≥5 methylated molecules |

*Median turnaround time from sample receipt to results.

[†]Limit of detection (LoD) defined as the allele fraction/copy number at which the test has a 95% probability of detection for oncogenic variants and genes with relevance in guidelines, drug labels, and clinical trials.

[‡]Indicates mutant allele fraction for detected SNVs and Indels.

[§]Tumor mutation burden (TMB) evaluable at or above a sample allele fraction of 0.3%.

CHIP: Clonal Hematopoiesis of Indeterminate Potential; CNA: Copy Number Amplification; MSI: Microsatellite Instability; SNV: Single Nucleotide Variant.

Important Note: The Guardant360 test was developed as a Laboratory Developed Test (LDT), and its performance characteristics determined, by the Guardant Health Clinical Laboratory in Redwood City, CA, USA, which is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) as qualified to perform high-complexity testing. This test has not been cleared or approved by the US FDA.



Assay Specifications

Gene Panel

MSI status - Qualitative result

TMB - Mutations per megabase

CHIP - Potential variants noted

Tumor fraction - Percent

| | | | | | | | | | | |
|-----------------------|------------------------|---------------------------|---------------------|----------|-----------------------|-----------------------|----------------------|------------------------|-----------------------|--------------------|
| ABCB1 | BTG1 | CUX1 | ERCC6L2 | FZD10 | INTS6L | MKNK1 | PBRM1 ^{††} | PTPRT | SHLD2 | TNK2 |
| ABL1 | BTG2 | CWC22 | ERCC8 | FZD2 | IRF1 | MLH1 ^{††} | PCBP1 | QKI | SLC34A2 | TNPO1 |
| ABL2 | BTK | CXCR4 | EREG | FZD3 | IRF2 | MLH3 | PCBP2 | RAB35 | SLFN11 | TOP1 |
| ABRAXAS1 | BUB1B | CYLD | ERF | FZD4 | IRF4 | MLST8 | PCDH15 | RAC1 | SLIT2 | TOP2A |
| ACVR1 | C9orf78 | CYP17A1 | ERG | FZD5 | IRS2 | MPL | PDCD1 | RAD18 | SMAD2 | TOPAZ1 |
| ACVR1B | CALR | CYP19A1 | ERRF1 | FZD6 | JAK1 | MRAS | PDCD1LG2 | RAD21 | SMAD3 | TP53 ^{††} |
| ADGRG4 | CARD11 | CYP2C19 | ESR1 ^{‡‡} | FZD7 | JAK2 | MRE11 | PDE7A | RAD50 ^{††} | SMAD4 ^{††} | TP53BP1 |
| ADARB2 | CASP8 | CYP3A4 | ETS1 | FZD8 | JAK3 | MSH2 ^{††} | PDGFRA [#] | RAD51 ^{††} | SMARCA2 | TP63 |
| ADGRA2 | CASR | DAXX | ETV1 | FZD9 | JUN | MSH3 | PDGFRB | RAD51B | SMARCA4 | TP73 |
| ADGRG4 | CAV1 | DCUN1D1 | ETV4 | GAS6 | KAT6A | MSH6 ^{††} | PDK1 | RAD51C ^{††††} | SMARCAL1 | TPMT |
| AFDN | CBFB | DDIT3 | ETV5 | GATA1 | KAT6B | MTAP ^{††††} | PDPK1 | RAD51D ^{††††} | SMARCB1 | TRAF2 |
| AGGF1 | CBL | DDR1 | ETV6 | GATA2 | KDM4A | MTHFR | PHF6 | RAD52 | SMARCD1 | TRAF3 |
| AIP | CBLB | DDR2 | EWSR1 | GATA3 | KDM5A | MTOR | PHLPP1 | RAD54L | SMARCE1 | TRAF7 |
| AKT1 | CCAR1 | DDX17 | EXO1 | GATA4 | KDM5B | MUTYH ^{††} | PHLPP2 | RAE1E | SMC1A | TRIM24 |
| AKT1S1 | CCN6 | DDX18 | EZH1 | GATA6 | KDM5C | MYB | PHOX2B | RAF1 [#] | SMC3 | TRIP13 |
| AKT2 | CCNA2 | DDX27 | EZH2 | GATA3 | KDM6A | MYC [#] | PIAS4 | RARA | SMO | TSC1 |
| AKT3 | CCNB1 | DDX3X | FAAP100 | GID4 | KDR | MYCL | PIK3C2B | RASA1 | SNCAIP | TSC2 ^{††} |
| ALB | CCND1 [#] | DDX41 | FAAP20 | GLI1 | KEAP1 ^{††††} | MYCN | PIK3CA [#] | RB1 ^{††††} | SOC1 | TSHR |
| ALK [#] | CCND2 [#] | DEPDC5 | FAAP24 | GNA11 | KIN | MYD88 | PIK3CB | RBBP6 | SOC3 | TSHZ2 |
| ALOX12B | CCND3 | DEPTOR | FANCA ^{††} | GNA13 | KIT [#] | MYO1D | PIK3CD | RBM10 | SOS1 | TYMP |
| ALOX15B | CCNE1 [#] | DHX15 | FANCB | GNAQ | KLF4 | NAB2 | PIK3CG | RBX1 | SOX10 | TYMS |
| ALOX5 | CCNE2 | DHX16 | FANCC | GNAS | KLHL6 | NBN | PIK3R1 | RECQL | SOX17 | TYRO3 |
| AMER1 | CD274 | DHX36 | FANCD2 | GPATCH8 | KLHL9 | NCOR1 | PIK3R2 | RECQL4 | SOX2 | U2AF1 |
| APC ^{††} | CD276 | DHX9 | FANCE ^{††} | GPC3 | KMT2A | NCR1 | PIK3R3 | RET [†] | SOX9 | UBE2T |
| APEX1 | CD74 | DICER1 | FANCF ^{††} | GREM1 | KMT2B | NCR3 | PIM1 | REV3L | SPEN | UGT1A1 |
| APLNLR | CD79A | DIS3L2 | FANCG ^{††} | GRIN2A | KMT2C | NEGR1 | PIN1 | RGS1 | SPOP | UIMC1 |
| AR [#] | CD79B | DLL4 | FANCI | GSK3B | KMT2D | NELFE | PKM | RHEB | SRC | ULBP1 |
| ARAF | CDC27 | DNAJB1 | FANCL ^{††} | GSTM1 | KNSTRN | NF1 ^{††} | PLCG2 | RHOA | SRSF2 | ULBP3 |
| ARFRP1 | CDC5L | DNMT1 | FANCM | GSTP1 | KRAS [#] | NF2 ^{††} | PLEKHS1 | RHOB | SRY | USP28 |
| ARHGAP35 | CDC7 | DNMT3A | FAS | H3-4 | LATS1 | NFE2L2 | PLRG1 | RICTOR | SS18 | USP7 |
| ARID1A ^{††} | CDC73 | DNMT3B | FAT1 | H3F3A | LGR4 | NFKBIA | PMS1 | RIF1 | STAG2 | USP9X |
| ARID1B ^{††} | CDH1 ^{††} | DOT1L | FBXW7 | HACD4 | LGR5 | NHEJ1 | PMS2 ^{††} | RILPL1 | STAT1 | VEGFA |
| ARID2 | CDH6 | DPYD | FCGR2A | HDAC2 | LGR6 | NKX2-1 | POLA1 | RIT1 | STAT3 | VEGFB |
| ASXL1 | CDK11A | DUSP4 | FCGR3A | HDAC6 | LIG1 | NOTCH1 | POLD1 | RNAHEH2B | STAT4 | VHL |
| ATM ^{††††} | CDK12 ^{††††} | DYNLL1 | FEN1 | HELQ | LIG4 | NOTCH2 | POLE ^{††} | RNF43 | STK11 ^{††††} | VIRMA |
| ATMIN | CDK4 [#] | DYRK2 | FGF1 | HES1 | LMO1 | NOTCH3 | POLH | ROBO1 | STK19 | WBP11 |
| ATR ^{††} | CDK6 [#] | E2F3 | FGF10 | HEY1 | LRP1B | NOTCH4 | POLQ | ROBO2 | STK40 | WEE1 |
| ATRX | CDK7 | ECT2L | FGF12 | HEYL | LRP2 | NOVA1 | POT1 | ROS1 [†] | STN1 | WRN |
| AURKA | CDK8 | EFTUD2 | FGF14 | HGF | LRP5 | NPM1 | POU2F2 | RPA1 | SUFU | WT1 |
| AURKB | CDKN1A | EGFR ^{‡‡} | FGF19 | HNF1A | LRP6 | NPRL2 | PPARG | RPS27A | SYK | WWP1 |
| AURKC | CDKN1B ^{††} | EIF1AX | FGF2 | HNRNPDL | LTK | NPRL3 | PIIG | RPS6KA3 | SYNCRIP | XBP1 |
| AXIN1 | CDKN1C | EIF4A1 | FGF23 | HOXB13 | LYN | NRAS | PPM1D | RPS6KB1 | TACSTD2 (TROP2) | XPA |
| AXIN2 | CDKN2A ^{††††} | EIF4A2 | FGF3 | HRAS | LZTR1 | NRG1 [†] | PPP2CA | RPS6KB2 | TAF1L | XPC |
| AXL | CDKN2B | EIF4A3 | FGF4 | HSD3B1 | MAD2L2 | NSD1 | PPP2R1A | RPTOR | TAP1 | XPO1 |
| B2M ^{††} | CDKN2C | EIF4B | FGF5 | HSP90AA1 | MALT1 | NSD2 | PPP2R2A | RRAGC | TAP2 | XRCC1 |
| BABAM1 | CEBPA | EIF4E | FGF6 | ICOSLG | MAP2K1 | NSD3 | PPP3CA | RSP01 | TABBP | XRCC2 |
| BABAM2 | CELF4 | EIF4E2 | FGF7 | ID3 | MAP2K2 | NSRP1 | PPP6C | RSP02 | TBC1D7 | XRCC3 |
| BAP1 | CEP295 | ELAVL1 | FGF8 | IDH1 | MAP2K4 | NTHL1 | PRDM1 | RSP04 | TBX3 | XRCC4 |
| BARD1 | CFAP20 | ELAVL2 | FGF9 | IDH2 | MAP3K1 | NTRK1 [†] | PREX1 | RUNX1 | TCERG1 | XRCC5 |
| BCL2 | CHD4 | ELF3 | FGFR1 ^{‡‡} | IDO1 | MAP3K13 | NTRK2 [†] | PREX2 | RUNX1T1 | TCF7L2 | XRCC6 |
| BCL2L1 | CHEK1 ^{††} | ELOC | FGFR2 ^{‡‡} | IFNG | MAP4K3 | NTRK3 [†] | PRKAR1A | RXRA | TEK | YAP1 |
| BCL2L2 | CHEK2 ^{††††} | EML4 | FGFR3 [†] | IFNGR1 | MAPK1 | NUMA1 | PRKCI | RYBP | TEN1 | YES1 |
| BCL6 | CIC | EMSY | FGFR4 | IFNGR2 | MAPK3 | NUMB | PRKDC | SAMHD1 | TENT5C | ZC3H13 |
| BCOR | CMTM4 | EP300 | FH | IFNW1 | MAPKAP1 | NUP93 | PRKN | SDC4 | TERT [†] | ZC3H18 |
| BCORL1 | CMTM6 | EPCAM | FLCN | IGF1 | MARK2 | NUTM1 | PRMT5 | SDHA ^{††} | TET1 | ZC3H4 |
| BCR | CNOT3 | EPHA3 | FLT1 | IGF1R | MAX | P2RY8 | PRPF40B | SDHAF2 | TET2 | ZMYM3 |
| BIRC5 | CREBBP | EPHA5 | FLT3 | IGF2 | MCL1 | PABPC1 | PRPF4B | SDHB ^{††} | TFE3 | ZNF217 |
| BLM | CRKL | EPHA7 | FLT4 | IGF2BP3 | MDC1 | PAK1 | PSENEN | SDHC ^{††} | TFRC | ZNF703 |
| BMPR1A | CRTC1 | EPHB1 | FOX A1 | IGF2R | MDM2 | PAK3 | PSMB10 | SDHD ^{††} | TGFBF1 | ZNF73 |
| BRAF ^{‡‡} | CSF1R | ERBB2 [#] (HER2) | FOX L2 | IKBKE | MDM4 | PALB2 ^{††††} | PSMB8 | SEM1 | TGFBF2 | ZNRF3 |
| BRCA1 ^{††††} | CSF3R | ERBB3 | FOX O1 | IKZF1 | MED12 | PARG | PSMB9 | SERPINB3 | THRAP3 | ZRSR2 |
| BRCA2 ^{††††} | CTC1 | ERBB4 | FOX P1 | IL1R1 | MEF2B | PARP1 | PTCH1 | SERPINB4 | TIA1 | |
| BRCC3 | CTCF | ERCC1 | FRS2 | IL2RA | MEN1 | PARP2 | PTDSS1 | SES2 | TIPARP | |
| BRD2 | CTLA4 | ERCC2 | FUBP1 | IL2RB | MERTK | PAX3 | PTEN ^{††††} | SETD2 | TMEM127 | |
| BRD3 | CTNNA1 | ERCC3 | FUBP3 | IL2RG | MET ^{††} | PAX5 | PTPN11 | SF3B1 | TMPPRSS2 | |
| BRD4 | CTNNB1 | ERCC4 | FUS | IL7R | MGA | PAX7 | PTPN2 | SF3B3 | TNFAIP3 | |
| BRIP1 ^{††} | CUL3 | ERCC5 | FYN | INHBA | MGMT ^{‡‡} | PAX8 | PTPRD | SH2D1A | TNFRSF14 | |
| BSG | CUL4A | ERCC6 | FZD1 | INPP4B | MITF | PAXIP1 | PTPRS | SHLD1 | TNFRSF1A | |

[†]Includes TERT promoter region. [‡]Includes CNAs. ^{‡‡}Includes Fusions/Rearrangements. ^{††}Includes Copy Number Losses. ^{†††}Includes Promoter Methylation. ^{††††}Promoter Methylation only. CHIP filter is based on a proprietary algorithm that identifies variants that are potentially derived from clonal hematopoiesis.